



FuseBox™ Processor

fusebox™

Video Inputs	DVI, Component (YPbPr), Y/C, Composite
Video Outputs	Digital: 3 x DVI, 3 x Component (YPbPr) Analog: 3 x Y/C, 3 x Composite (NTSC)
Communication Interfaces	USB, RS-232, Touch panel
Power	Voltage: 115 VAC / 230 VAC Frequency: 50Hz / 60Hz Voltage Fluctuation: +/- 10% Max Power: 300W
Insufflation Interface	Irrigation and optional CO ₂
Dimensions (W x H x D)	350mm x 170mm x 450mm (13.78" x 6.7" x 17.72")
Weight	9 kg (19.8 lb)
Electrical Class	Class I
Electrical Type	Type BF
Air Pump Type	Diaphragm Max: 0.45 Bar Irrigation flow rate: 1.8cc/Sec

Includes

Leak Tester

Water Bottle (2)

White Balance Cap