



FuseCart™ Endoscopy Carts

**fuse**cart™



## Deluxe Motorized FuseCart™

### EDM-100P

#### Dimensions (W x H x D)

1.07m x 1.08mm x .6m  
(42" x 42.5" x 23")

Designed by Endoscopy Nurses as the most EFFICIENT and SAFE travel system

Safety features include variable speed, automatic braking, directional control, and dual handles

Ideal working and driving height, optimizing your Endoscopy Workstation

Retractable center drive wheel allows easy maneuvering in tight spaces

Customizable with a variety of storage features (baskets, drawers, shelving, scope hanger, etc)

Flexible monitor mounting including single, dual, and docking monitor stand

Comes standard with the following: 2 shelves, 3 drawers, fold-down keyboard tray and mouse pad, and power strip

## Deluxe FuseCart™

### ED-100P

#### Dimensions (W x H x D)

1.07m x 1.08mm x .6m  
(42" x 42.5" x 23")

Designed by Endoscopy Nurses to optimize your Endoscopy Workstation

Safety features include two directional steering casters, dual handles, and two braking casters

Ideal working height for the procedure room or mobile system

Reinforced aluminum frame provides lightweight durability

Flexible monitor mounting including single, dual, and docking monitor stand

Customizable with a variety of storage features (baskets, drawers, shelving, scope hanger, etc)

Comes standard with the following: 2 shelves, 3 drawers, fold-down keyboard tray and mouse pad, and power strip

## Standard FuseCart™

### ES-100P

#### Dimensions (W x H x D)

.6m x 1.08mm x .6m  
(22.5" x 42.5" x 23")

Compact and efficient with raised-edge working surface, ideal for any procedure room

Customizable with a variety of storage features (baskets, drawers, shelving, scope hanger, etc)

Ideal working height for the procedure room

Reinforced aluminum frame provides lightweight durability

Flexible monitor mounting including single or dual flat panel monitors

Comes standard with: 2 shelves, fold-down keyboard tray with mouse pad, and power strip

