



fuse

ENDOCHOICE

Fuse 1C™ Colonoscope

fuse

Insertion Tube OD	12.8mm
Distal Tip OD	13.9mm
Working Channel	3.8mm
Working Length	1680mm
Field of View	330°
Depth of Field	100mm
Tip Deflection Up/Down	180°/180°
Tip Deflection Right/Left	160°/160°
Control Buttons	Capture, Zoom, Record, Freeze, Timer
Illumination Method	LED

Includes

Endoscope Case

Endoscope Reprocessing Adaptor Kit

Reprocessing Manual