



**Fuse 1G™ Gastroscope**

**fuse™**

<b>Insertion Tube OD</b>	9.6mm
<b>Distal Tip OD</b>	10.5mm
<b>Working Channel</b>	2.8mm
<b>Working Length</b>	1050mm
<b>Field of View</b>	220°
<b>Depth of Field</b>	100mm
<b>Tip Deflection Up/Down</b>	210°/120°
<b>Tip Deflection Right/Left</b>	120°/120°
<b>Control Buttons</b>	Capture, Zoom, Record, Freeze, Timer
<b>Illumination Method</b>	LED

### Includes

Endoscope Case

Endoscope Reprocessing Adaptor Kit

Reprocessing Manual